

HEALING IS VOLTAGE™ **DAILY PROTOCOL**

Please make sure you do the steps in order.





TABLE OF CONTENTS

Table of Contents	1
Step 1: Recharge Primary Voltage Circuit	2
Step 2: Reset Autonomic Nervous System	3-4
Bonus: Spinal Reset	5
Step 3: Correct Polarity	6
Step 4: Activate Fetal Wiring System	7
Step 5: Full Body Recharge	8-12
Bonus: Full Body Charge with Pads	13
Treat Areas of Concern	14



Scan this QR code to watch a video walk-through of the Tennant Daily Protocol.

**Questions? Schedule an appointment
with a Senenergy Advisor today.**
www.senenergy.us/appointments

STEP 1

RECHARGE THE PRIMARY VOLTAGE CIRCUIT (LIVER/GALL BLADDER MERIDIAN)



Place the Y-Electrode right at the base of the hairline. This can be done on your own, or with someone holding the Y-Electrode in place.

SETTINGS:

- Connect BioModulator® to Y-Electrode (or, use electrode on back of BioModulator®).
- Set BioModulator® to Ten-8 Mode.
- Hold electrode at the base of the hairline on the back of the neck and turn power up until you feel a tingle.
- **PRO:** Hold until the display reads DZ.
- **PLUS:** Hold for 3 minutes.



Why is this step important?


The cranial sacral pump that runs from our head to the base of our spine transports cerebral spinal fluid back and forth to carry waste out of the brain and nutrients in.

When we apply voltage to this location, it not only helps us switch our nervous system into a parasympathetic (rest & digest) state and ensure everything is flowing properly, but it also helps close some of the digestive valves that stay open when stuck in a sympathetic (fight or flight) state, preparing your body for the next steps.

STEP 2

RESET AUTONOMIC NERVOUS SYSTEM (BOWLING BALL SYNDROME)

*For managing anxiety, depression,
insomnia, stress, and PTSD.*



STEP 2

This should be relaxing, like a mini meditation. Find a quiet place to do this therapy. Be patient, close your eyes and relax.

Place the Y-Electrode (or BioModulator® back electrode) on the muscle below the ear. Make sure that you have the electrode firmly pressed against the SCM muscle as pictured to the right.

The X marks the back of the muscle where the electrode is placed as shown in the picture to the right and on the previous page.

Select Ten-8 mode on the BioModulator® and turn up the power until you feel a gentle tingling sensation.


Note: You will select the power level in Ten-8 mode (a constant frequency) to ensure your comfort, but you will be doing the therapy in Infinity, which is a modulation frequency. You will feel the frequency intermittently, which is why you don't want to set your comfort level on Infinity.

Again, once you select your comfortable power level in Ten-8, **switch** to Infinity mode **without taking the electrode off the neck.**

Hold for 1-3 minutes per side.

Please note: You will know you are receiving the therapy when you feel a strange sensation going up the ear or jaw.

This can be done as many times a day as necessary. If you do this step right before bed, your body will get better rest!



SCM muscle placement

Place the Y-Electrode (or BioModulator® back electrode) on the muscle below the ear.

Why is this step important?

By stimulating the vagus nerve with voltage, you signal the body to turn the parasympathetic nervous system back on, restoring full function to the organs of your body and reducing overall stress.

This is an additional bonus step if you have a partner:

Reset the autonomic ganglia along the spine using Ten-8 Mode on the BioModulator® PLUS or PRO.

Place the Y-Electrode at the top of your spine, set the BioModulator® to **Ten-8 Mode** and turn the power up until you gently feel it.

Move the electrode **down the spine slowly from top to bottom 10 times.**

(Do not move from bottom to top)



**This corrects the
polarity of all circuits.**

SETTINGS:

Direct the BioTransducer® either up the spine between the legs or hold it above your head with the device pointing down the spine, to correct the circuits and BioTerminal points.




For Pro, set mode to
783 Harmonics
and set power to **250**.

OR



For Plus, set mode to
Ten-8 and set power to
the **solid orange dot**.



Duration: **5 minutes**
for **BioTransducer® Pro**



OR

Duration: **15 minutes**
for **Crystal Wave**



STEP 3

CORRECT POLARITY



If you don't have a stand to assist you in holding it on top of your head, you may complete this step by lying down and resting the BioTransducer® on the crown of the head or between the legs. If it is on top of your head, shine the BioTransducer down your spine, and if it is between your legs, shine it up towards your spine.

Why is this step important?

By ensuring our body's polarity is pointing in the right direction, this allows our muscle battery packs to hold a charge when we apply voltage in Step 5.

STEP 4

ACTIVATE FETAL WIRING SYSTEM



Point the BioTransducer® at the belly button. This step can be done with clothes on. The energy emitted from the BioTransducer® will pass through clothing.

This will turn on all the organs of the body's original wiring system found at the umbilicus.

SETTINGS:

Connect the BioModulator® to the BioTransducer® and set to **Infinity mode**, then raise the power of the BM to the **highest level**.



PRO: Set mode to **Infinity** and turn power up to **250**.

OR



PLUS: Set mode to **Infinity** and turn power up to the **solid orange dot**.



DURATION:



PRO: 5 minutes per day


OR



Crystal Wave: 15 minutes per day

STEP 5

TOTAL BODY RECHARGE



You may perform this step sitting up or laying down. Also, it is normal to only feel the tingling sensation on one or both hands during this step of the therapy. The soles of the feet are thicker, so you may not feel any tingling there.

This recharges the whole body's battery packs. Remember you have to correct the polarities before the batteries will take a charge!

If you are active/using for preventative you can skip this step.

SETTINGS:

Connect the handgrips and footplates to the BioModulator®.



PRO:
Solfeggio Mode



PLUS:
Ten-8 Mode

OR

Left hand **RED** Right hand **BLACK** Left foot **BLACK** Right foot **RED**

While holding the hand grips, power the BioModulator® to what is **gentle and comfortable**.

Duration: **30 min**

See next pages for how to configure the hand grips and foot plates.

HAND GRIP + FOOT PLATE CONFIGURATION

STEP 1

Grab the separate **red** and **black** banana bridge wires.

These will “bridge” the hand grips and foot plates together so you can use them at the same time.



STEP 2

Connect one end of the **red** bridge wire to the **red** foot plate.

-and-

Connect one end of the **black** bridge wire to the **black** foot plate.

CONFIGURATION CONTINUED...

STEP 3

Connect the other end of the **red** bridge wire to the **red** hand grip.

-and-

Connect the other end of the **black** bridge wire to the **black** hand grip.



STEP 4

Grab the lead wire that connects to the BioModulator® and has an end with a **red** and **black** split wire.


Connect to the BioModulator®.



CONFIGURATION CONTINUED...


STEP 5

Insert the **red** end of the lead wire into the open space on the **red** bridge wire up near the hand grip.



STEP 6

Insert the **black** end of the lead wire into the open space on the **black** bridge wire up near the handgrip.



CONFIGURATION CONTINUED...

STEP 7

Hand grips and foot plates are now configured to work together at the same time.



REMEMBER:

LEFT HAND
RED

RIGHT HAND
BLACK

LEFT FOOT
BLACK


RIGHT FOOT
RED

ALL SET

TOTAL BODY RECHARGE (with pads)

FRONT OF BODY

BACK OF BODY



This is a great option if you didn't travel with your hand grips & foot plates, or if you'd like to multitask while charging.

SETTINGS:



For the BioModulator® Plus & Pro, use **Ten-8 Mode.**



Please apply the conductive pads (patches) where you see the red and black dots illustrated on the front & back of the body (four pads total).

Necessary wires/pads:



**4PIN Connection 2/2
Lead Wire**

(4 pads at a time)

Red/Black Wire = Red/Black Patch



3" Round Foam Pads

You must wear the patches for at least 4-6 hours. You may wear them longer as desired. For best results, the pad area must be clean and bare (no hair).



TREAT AREAS OF CONCERN

Now that you have completed the daily steps, your body is now primed and prepped to receive additional therapies on areas of concern like general pain, skin therapy, organ therapy, etc.

Use the BioModulator®'s built-in electrode, or connect it to your Y-Electrode or conductive pads of choice.

Select Ten-8 Mode to treat pain in the muscles and bones, and Infinity Mode to treat organs.

If you have a BioModulator® Pro, you may use one of the 60 modes that is better suited to your area of concern.

Feel free to experiment and switch between the modes to see which one your body responds to the most.

Using the BioTransducer®:

- Connect your BioTransducer® to your BioModulator®
- BioModulator® PRO: Increase power to 250
- BioModulator® PLUS: Increase power to the solid orange dot
- Point the BioTransducer® at the area for 15-20 min of therapy.

It is important to listen to your body and make modifications if needed.